

HOW TO PUT ON, USE, TAKE OFF AND DISPOSE OF A MASK

HOW TO PUT ON, USE, TAKE OFF AND DISPOSE OF A MASK

Before putting on a mask,
clean hands with alcohol-based hand rub or soap and water



HOW TO PUT ON, USE, TAKE OFF AND DISPOSE OF A MASK

Cover mouth and nose
with mask and **make sure there are no gaps between your face and the mask**




HOW TO PUT ON, USE, TAKE OFF AND DISPOSE OF A MASK

Avoid touching the mask while using it;
if you do, clean your hands with alcohol-based hand rub or soap and water



HOW TO PUT ON, USE, TAKE OFF AND DISPOSE OF A MASK

Replace the mask with a new one as soon as it is damp and do not re-use single-use masks



HOW TO PUT ON, USE, TAKE OFF AND DISPOSE OF A MASK

To remove the mask: **remove it from behind** (do not touch the front of mask); **discard immediately in a closed bin;** **clean hands with alcohol-based hand rub or soap and water**

